

Converting Volumetric Flow Rate (ml/min) to Linear Flow Rate (cm/hour)

$$\text{Linear flow rate (ml/min)} = \frac{\text{Volumetric flow rate (ml/min)} \times 60}{\text{Cross sectional area of column (cm}^2\text{)}}$$

$$= Z \times 60 \times \frac{4}{\pi \times d^2}$$

Z = volumetric flow rate in ml/min

d = column inner diameter in cm

Table for Converting Volumetric Flow Rates (ml/min) to Linear Flow Rates (cm/h) for Columns of Varying Internal Diameters (mm):

Choose the volumetric flow rate that you are using against the internal diameter of your column to find the equivalent flow rate in cm/h.

For example a WorkBeads™ 40 IEX 8mm pre-packed column (8mm internal diameter) will have a Linear Flow Rate of 238.85 cm/h at 2.0 ml/min.

		Internal Diameter of Column (mm)								
		5.0	8.0	10.0	12.5	16.0	26.0	50.0	100.0	180.0
Vol. Flow Rate (ml/min)	0.1	30.57	11.94	7.64	4.89	2.99	1.13	0.31	0.08	0.02
	0.2	61.15	23.89	15.29	9.78	5.97	2.26	0.61	0.15	0.05
	0.3	91.72	35.83	22.93	14.68	8.96	3.39	0.92	0.23	0.07
	0.5	152.87	59.71	38.22	24.46	14.93	5.65	1.53	0.38	0.12
	1	305.73	119.43	76.43	48.92	29.86	11.31	3.06	0.76	0.24
	1.5	458.60	179.14	114.65	73.38	44.79	16.96	4.59	1.15	0.35
	2	611.46	238.85	152.87	97.83	59.71	22.61	6.11	1.53	0.47
	3	917.20	358.28	229.30	146.75	89.57	33.92	9.17	2.29	0.71
	4	1222.93	477.71	305.73	195.67	119.43	45.23	12.23	3.06	0.94
	5	1528.66	597.13	382.17	244.59	149.28	56.53	15.29	3.82	1.18
	10		1194.27	764.33	489.17	298.57	113.07	30.57	7.64	2.36
	20			1528.66	978.34	597.13	226.13	61.15	15.29	4.72
	30				1467.52	895.70	339.20	91.72	22.93	7.08
	40					1194.27	452.27	122.29	30.57	9.44
	50					1492.83	565.33	152.87	38.22	11.80
	100						1130.67	305.73	76.43	23.59
	200							611.46	152.87	47.18
300							917.20	229.30	70.77	
500								382.17	117.95	